

High-performance analysis of petrochemicals by GC×GC with parallel detection by FID and MS

This study demonstrates the performance of a flow-modulated GC×GC–MS system for the analysis of complex petrochemicals, with the INSIGHT modulator providing efficient and repeatable modulation without the need for expensive liquid cryogen. Analyte speciation is also aided by Tandem Ionisation for simultaneous acquisition of hard and soft EI data.



Introduction

Precise characterisation of petrochemical samples is crucial for quality control, and also to understand the reactions that take place during refining processes. Comprehensive 2D gas chromatography (GC×GC) offers significant advantages over 1D chromatography for such analyses, with its vastly expanded separation space and the added benefit of highly structured groupings of compounds.

The key component in the GC×GC system is the modulator – the device that samples and re-injects the first-column effluent on to the second column in narrow bands to ensure that the first-dimension separation is retained and that the short second-dimension column does not become overloaded.

Thermal modulation is the most commonly used technique, but this often requires expensive liquid cryogen and can make it difficult to achieve precise replication of results across multiple instruments. An alternative approach is flow modulation, which avoids the inconvenience and expense of liquid cryogen, and offers much better between-sample and between-instrument repeatability.

Fast acquisition speeds and sensitive detection make time-of-flight mass spectrometry (TOF MS) the ideal partner to GC×GC. However, best-practice in the petrochemical industry recommends the use of flame ionisation detection (FID) for target compound quantitation. This study explores the use of parallel-detection GC×GC to enable both targeted and non-targeted approaches to be performed simultaneously – with TOF MS for confident identification and FID for robust quantitation.

Experimental

A schematic of the analytical system is shown in Figure 1.

Sample: Pump diesel (undiluted).

GC×GC: Injector:

Split/splitless; Liner:

Single taper with

wool, 4 mm (i.d.); Carrier gas: Helium, constant-flow at 0.5 mL/min; Mode: Split

300:1; Temperature: 280°C; Septum purge: On, 3 mL/min. Modulator: INSIGHT

flow modulator (SepSolve Analytical). **Repeatability study:** 2D column set: 1st

dimension: BPX5, 20 m × 0.18 mm × 0.18 μm; 2nd dimension: DB17, 2 m ×

0.25 mm × 0.15 μm. Temperature program: Main oven: 50°C (1 min), 3°C/min to

320°C (hold 10 min). Loop dimensions: 30 cm × 0.32 mm i.d.; Loop volume:

24 μL; Fill time: 5800 ms; Flush time: 200 ms; P_M : 6.0 s. **Dual detection study:** 2D

column set: 1st dimension: BPX5, 25 m × 0.15 mm × 0.25 μm; 2nd dimension:

ZB-50, 5 m × 0.25 mm × 0.15 μm. Temperature program: Main oven: 50°C

(1 min), 2.5°C/min to 320°C (hold 20 min). Loop dimensions: 10 cm × 0.53 mm i.d.

(loop volume: 22 μL); Fill time: 3900 ms; Flush time: 100 ms; P_M : 4.0 s.

TOF MS: Instrument: BenchTOF-Select™ (Markes International); Filament voltage: 1.7 V; Ion source: 320°C; Transfer line: 300°C; Mass range: m/z 45–450; Data rate: 100 Hz in Tandem Ionisation® mode at 70 eV and 14 eV.

FID: H₂ flow: 30 mL/min; N₂ flow: 25 mL/min; Air flow: 400 mL/min; Data rate: 100 Hz.

Software: Instrument control and GC×GC data processing was carried out using TOF-DS™ with ChromSpace®.




Figure 1

Schematic of the GC×GC-TOF MS/FID analytical system.

Results and discussion

Repeatability of flow modulation

A major advantage of flow-modulated GC×GC over thermally-modulated systems (aside from the reduction in running costs) is the superior repeatability that can be achieved. INSIGHT's precisely-defined microfluidic design^[1] allows identical configurations to be installed across multiple instruments, unlike in thermal devices, where small variations in column position can have a large impact on results.

In this study, the repeatability of the flow modulation device was first evaluated using the GC×GC-FID results for 24 replicate injections over the course of three days (Figure 2).

The charts in Figure 2 illustrate the excellent repeatability of the technique, with calculated RSDs less than 0.14%, 0.52% and 2.98% for 1t_R , 2t_R and peak abundance, respectively. This high precision of time response is aided by having a dedicated EPC unit for each column, and enables simple comparison of large sample batches and confident interpretation of results.

Evaluation of GC×GC with parallel MS/FID detection

FID is the gold standard for quantitative GC and GC×GC applications in the petrochemical industry for many reasons, including negligible band broadening, wide linear range, fast data capture and mass-dependent response. At the same time, MS detection (especially high-performance TOF MS), offers the advantage of robust analyte identification, and options for peak deconvolution and retrospective data-searching.

The use of GC×GC with parallel detection by MS and FID is therefore an ideal approach for achieving qualitative and quantitative analysis of petrochemicals in a single run. Targeted quantitative analysis can be performed across the widest concentration ranges, with dual-detector cross-validation and confirmation of peak purity.




Figure 2

Repeatability of GC×GC-FID using the INSIGHT modulator during a three-day study (n = 8 per day). The RSDs for each day are indicated.

However, thermally-modulated systems can be troublesome to configure for such dual-detector setups, but this is not an issue with flow modulation. Figure 3 compares colour plots obtained by simultaneous detection by TOF MS and FID, using chromatographic conditions for optimum separation and alignment of retention times in both the first and second dimensions.




Figure 3

GCxGC colour plots of diesel, generated by simultaneous detection by TOF MS (with Tandem Ionisation at 70 eV and 14 eV) and FID.

Figure 4 further illustrates the excellent retention time alignment, using the linear TOF MS and FID traces.




Figure 4


Overlays of the FID (blue) and TOF MS (black) linear traces, showing the excellent retention time alignment achieved throughout the run, using the INSIGHT modulator.

Tandem Ionisation^[2]

The use of simultaneous detection by TOF MS and FID provides a wealth of information for petrochemical analyses, but the challenge of identifying individual hydrocarbon isomers still remains. At conventional 70 eV ionisation, the spectra of these isomers are virtually indistinguishable, meaning it is impossible to make confident compound assignments. Such identifications are of great commercial importance to the petrochemical industry, because individual isomers can provide different characteristics to the final fuel (e.g. octane number, flash point and viscosity).

Addressing this challenge, the use of soft electron ionisation (or soft EI) has been shown to provide increased confidence in the identification of branched alkanes,^[3] by providing enhancement of the molecular ion and other structurally-significant ions. A further advance in this area is Tandem Ionisation, which involves fast switching between conventional 70 eV ionisation and Select-eV[®] soft EI. This allows two complementary MS datasets to be simultaneously acquired, without impacting laboratory workflow and with perfectly aligned peaks for easy navigation of the soft EI data.

Figure 5 shows the speciation of three hydrocarbon isomers using soft EI. The 14 eV spectra exhibit clear spectral differences that enable confident identification – which would not have been achievable based solely on 70 eV analyses.

**Figure 5**

Spectral comparisons at 70 eV and 14 eV for three hydrocarbon isomers.

Conclusions

This study has illustrated the power of flow-modulated GC×GC to provide structured chromatograms that greatly simplify the process of compound identification, as well as the ability of the INSIGHT™ flow modulator to provide robust, repeatable analysis of petrochemicals, without the inconvenience or running costs associated with thermal modulation. In addition, the excellent retention-time repeatability in both dimensions enables fast and confident data processing of large sample batches, while the retention-time alignment between FID and TOF MS during parallel detection enables two-fold validation of measured peaks.

The analytical approach detailed here has also been shown to push the boundaries of hyphenated systems by providing three information-rich datasets per acquisition – namely, robust quantitation by FID, untargeted screening by TOF MS with electron ionisation (EI) at 70 eV, and confirmatory identification by TOF MS with soft EI.

For more information on this application, or any of the techniques or products used, please contact SepSolve.

References and notes

- [1] The valve-based INSIGHT modulator uses differential flows to 'fill' and 'flush' a sample loop (Figure 6) – meaning that running costs for routine GC×GC are low, and that drawbacks associated with liquid cryogen are avoided.




Figure 6

Flows used in the INSIGHT modulator to fill and flush the sample loop (orange).

- [2] Select-eV capability, available on the BenchTOF-Select mass spectrometer from Markes International, allows soft EI spectra to be collected down to 10 eV, and is fully automated by the instrument's software with no inherent loss in sensitivity or need for manual intervention. Tandem Ionisation now allows soft and hard ionisation spectra for a single peak to be simultaneously obtained, in both GC and GC×GC analyses, enabling challenging compounds (such as structurally similar isomers) to be discriminated without impacting laboratory workflows. Contact SepSolve for more details.
- [3] M.S. Alam, C. Stark and R.M. Harrison, Using variable ionisation energy time-of-flight mass spectrometry with comprehensive GC×GC to identify isomeric species, *Analytical Chemistry*, 2016, 88: 4211–4220, <http://dx.doi.org/10.1021/acs.analchem.5b03122>.

BenchTOF™, BenchTOF-Select™ and TOF-DS™ are trademarks of Markes International. ChromSpace®, Select-eV® and Tandem Ionisation® are registered trademarks of Markes International.

INSIGHT™ is a trademark of SepSolve Analytical.

Applications were performed under the stated analytical conditions. Operation under different conditions, or with incompatible sample matrices, may impact the performance shown.

WP006_2_071016