

GERSTEL Solutions

Completely Automated Hydrolysis, Extraction and Analysis
of Opioids in Urine using a New Robotic Autosampler and
LC/MS/MS Platform

Fred Foster


Applications Scientist

GERSTEL, Inc.

Baltimore, MD, USA




3 Waves of the Rise in Opioid Overdose Deaths



SOURCE: National Vital Statistics System Mortality File.

Potency Compared to Morphine




<https://www.washingtonpost.com/graphics/2017/health/opioids-scale>

Opioid Testing

Matrix


- Blood
- Plasma
- Serum
- Urine
- Dried Blood Spots
- Saliva




Analytical Technique

ELISA, GC/MS, LC/MS/MS


Morphine




Morphine-6-B-D-glucuronide




Hydromorphone




Hydromorphone-3-B-D-glucuronide




Codeine




Codeine-6B-glucuronide



Buprenorphine



Buprenorphine-3B-glucuronide



Hydrolysis – acid hydrolysis or enzymatic hydrolysis

“Degradation of Opioids and Opiates During Acid Hydrolysis Leads to Reduced Recovery Compared to Enzymatic Hydrolysis”


- P. Sitasuwan, et al, Journal of Analytical Toxicology, Volume 40, Issue 8, 1 October 2016, Pages 601–607

Variety of B-glucuronidase enzymes:

- A) Helix-Pomatia
- B) Abalone
- C) Recombinant Protein



Urine Samples treated with different β -Glucuronidase enzymes.





Extraction

- Dilute and shoot
- Precipitation (“solvent crash”)
- Liquid-Liquid extraction
- Solid Phase Extraction



Solid Phase Extraction





MPS

GERSTEL

MultiPurpose Sampler (MPS)

Integrated with LC/MS/MS



MPS

GERSTEL

MultiPurpose Sampler (MPS)

Benchtop Workstation



Integrated with GC/MS



Evaporative
Concentration **mVAP**



Extraction
Derivatization
Addition of Standards



Solid Phase
Extraction **SPE**
& Filtration



Cold storage
of samples



MPS for LC & LC/MS




CF200 Centrifugation



DisposablePipette-
Extraction **DPX**



Dried Blood Spot
Autosampler **DBS-A**



Online-SPE **SPE^{xos}**



Sample Preparation

Typical Manual Sample Preparation Techniques


- Pipetting
- Dilution
- Vortexing/Mixing
- Heating
- Sonication
- Evaporation
- Calibration Standard Preparation
- Liquid-Liquid Extraction
- Solid Phase Extraction


Automated Disposable Pipette Extraction DPX




DPX

GERSTEL DPX automation







LC Method Parameters

Mobile Phase:

A - 5mM ammonium formate with 0.05% formic acid
B –0.05% formic acid in methanol

LC Gradient:

Time (min)	Flow (mL/min)	Pressure (bar)	% B
0	0.5	800	5
0.5	0.5	800	5
1.5	0.5	800	30
3.5	0.5	800	70
4.5	0.5	800	95
6.49	0.5	800	95
6.5	0.5	800	5

Run time:

8 minutes.

Injection volume:

2.0uL (loop over-fill technique)

Analytical Column:

Agilent Poroshell 120 EC-C18 column, (3.0 x 50mm, 2.7um)

Column Temperature:

55 °C

Mass Spectrometer Parameters:

Operation:

Electrospray positive mode (Jet Stream)

Gas Temperature:

350°C

Gas Flow (N₂):

5 L/min

Nebulizer pressure:

35 psi

Sheath Gas Flow (N₂):

11 L/min

Sheath Gas Temperature:

250 °C

Capillary voltage:

4000 V

Nozzle voltage:

500 V



Delta EMV:

0 V


Compound Name	Precursor Ion (m/z)	Product Ion (m/z)		Frag (V)		CE (V)		Ret Time (min.)	High Std Conc. (ng/mL)	MRL (ng/mL)	Limit of Quantitation (ng/mL)
Buprenorphine ⁴	468.3	396.2	55.1	200	200	39	58	4.08	200	10.0	5.00
Codeine ³	300.2	165.1	128	158	158	45	60	2.33	1000	50.0	25.0
D3-morphine ¹	289.0	165.1	152	153	153	38	66	1.09	-	-	-
D3-Norbuprenorphine ²	417.3	152	55.1	190	190	124	76	3.62	-	-	-
D3-Tramadol ³	268.2	58.1	-	102	-	14	-	3.16	-	-	-
D4-Buprenorphine ⁴	472.3	400.2	59.1	210	210	40	56	4.06	-	-	-
D4-Meperidine ⁵	252.2	224.2	178.2	140	140	15	15	3.33	-	-	-
D5-Fentanyl ⁶	342.3	188.1	105.1	92	92	20	40	3.67	-	-	-
D5-Norfentanyl ⁷	238.2	84.1	55.1	125	125	20	50	3.06	-	-	-
D9-Methadone ⁸	319.3	268.1	-	118	-	8	-	4.31	-	-	-
Fentanyl ⁶	337.2	188.1	105.1	143	143	17	37	3.68	20.0	1.00	0.500
Furanyl Fentanyl ⁶	375.2	188.1	105	150	150	15	25	3.67	20.0	n/a	0.500
Hydrocodone ³	300.2	199	128	159	159	27	63	2.70	1000	50.0	25.0
Hydromorphone ¹	286.2	185	157	159	159	27	43	2.20	1000	50.0	25.0
Meperidine ⁵	248.2	220.1	174.1	128	128	17	13	3.33	1000	50.0	25.0
Methadone ⁸	310.2	265.1	105	112	112	7	27	4.32	1000	50.0	25.0
Morphine ¹	286.2	165.1	152	158	158	45	64	1.10	1000	50.0	25.0
Norbuprenorphine ²	414.3	187.1	83.1	205	205	37	53	3.63	200	10.0	5.00
Norfentanyl ⁷	233.2	84.1	55.1	112	112	12	36	3.08	20.0	1.00	0.500
Oxycodone ³	316.2	298.1	241.1	143	143	13	25	2.42	1000	50.0	25.0
Oxymorphone ¹	302.1	227.1	198.1	133	133	24	44	1.29	1000	50.0	25.0
Sufentanil ⁸	387.2	238.1	110.9	145	145	10	20	4.07	20.0	n/a	0.500
Tramadol ³	264.2	58.1	42.1	107	107	14	106	3.16	500	25.0	12.5

LC/MS Tool


GERSTEL





- Near-Zero carryover
- Low detection limits
- Fast
- Reproducible
- Integrated MAESTRO software control








GERSTEL



Resulting QC Accuracy and Precision



Drugs of Abuse in Oral Fluids: Automated SPE Extraction and LC/MS/MS Determination using a Robotic Autosampler

1. GERSTEL, Inc., 701 Digital Drive, Suite J, Linthicum, MD, 21090, USA
2. ITSP Solutions, Inc., 10 South Carolina St, Hartwell, GA 30643, USA

ABSTRACT Solid phase extraction (SPE) is one of the sample preparation methods most widely used by chromatographers, as can be seen from the huge number of SPE methods found in the literature. Typically, a liquid sample is passed across an adsorbent bed to retain and concentrate target analytes and eliminate interference from the sample matrix. Alternatively the adsorbent can be used to retain interferences while allowing the target analytes to pass through. Manual SPE cartridge formats can vary from disks to 96-well cartridges with various volumes of solvent to be tedious. However, solid phase extraction can be tedious and time consuming when performed manually and there is increasing demand for automation of SPE methods.

The MPS robotic[™] is a highly efficient LC or GC autosampler with extended robotic functionality. The MPS robotic[™] provides reliable processing of complex tasks including automation of SPE procedures. Syringe holders and syringes are integrated in special syringe modules, which can be exchanged automatically within a running sequence for maximum flexibility. The system is controlled in a simple and efficient manner using the proven GERSTEL MAESTRO software. In

this study, we show that an existing SPE method using ITSP cartridges [1] can be automated using the MPS robotic[®] autosampler under MAESTRO control. An example of a solid phase extraction method illustrating the determination of drugs of abuse in oral fluid is shown.

INTRODUCTION

INTRODUCTION
Drug testing based on the analysis of oral fluid samples is increasingly used to sample collection advantages over other biological matrices and due to the early detection window, it offers compared with other approaches. The collection of oral fluid is noninvasive and can be performed under directly observational conditions. The simplicity of collecting an oral fluid sample also reduces the likelihood of test sample adulteration or sample swapping leading to false positive or false negative results. Some collection devices, including the Quantalabs devices (2), used within this study, have an integrated seal to ensure an exact volume of oral fluid is collected. Following the collection step, the devices are subsequently placed in storage tubes containing stabilization buffer, ensuring that the sample can be shipped or stored for subsequent laboratory analysis.

While the buffer helps to stabilize analytes present in the oral fluid sample, it may also cause analyte suppression effects in the LC/MS/MS analysis leading

ter, Oscar Cabrices, Jackie Whitecavage,
ff, Edward Pfannkoch
o., 701 Digital Dr. Suite J,
MD 21090, USA

DS
reparation, Lab Automation, LC/MS/MS

Liquid extractions have long been performed manually used to extract and concentrate analytes from matrices. Inclusion of liquid-liquid extraction in analytical methods attests to the wide acceptance of the technique. Following solvent extraction it is also common to remove the solvent by evaporation and reconstitution step to improve sensitivity or exchange solvents for compatibility with chromatographic separations. Modern analytical methods are moving toward automation to help reduce solvent usage, increase sample throughput while ensuring the high quality of the resulting data.

Sampler (MPS), commonly used for sample collection in GC or HPLC can be used to perform a wide range of sample preparation techniques using a single autosampler and controlling software. The sampler can be interfaced to a GC or LC system or can be configured for stand-alone workstation and can also include a six position sample tray.

Thus, the automation of liquid-liquid extractions using an autosampler is discussed. Examination of a vortexing option that allows samples to be actively mixed using speeds of up to 3000 rpm compared liquid-liquid extractions methods for samples from different matrices are examined. Precision and accuracy data are provided.



Fred D. Foster, John R. Stuff, Edward A. Pfannkoch
GERSTEL, Inc., 701 Digital Dr. Suite J,
Linthicum, MD 21090, USA

KEYWORDS

KEYWORDS: Forensics, Veterinary, Sample Preparation, LC/MS/MS, High Throughput Lab Automation

ABSTRACT

The extraction of dried blood spots (DBS) typically involves manual intervention. First, a small disc is punched out of the center of a dried blood spot placed on a DBS card. Following solvent extraction of the sample, it is also common to include further cleanup steps, using solid phase extraction (SPE) to improve detection limits or exchanging solvents for compatibility with subsequent chromatographic separations. Modern analytical labs are looking to automate the process to help reduce solvent usage and to increase sample throughput while ensuring the high quality of the resulting data.

while ensuring the high quality and true results of the analysis. In this study, the use of dried blood spot analysis is demonstrated and the results evaluated. A novel, automated DBS Autosampler (DBS A) automatically inserts DBS cards into a Flow Through Desorption (FTD)[®] cell in which individual blood spots are rapidly and effectively desorbed. The DBS A is integrated into a complete clean-up and analysis system using online SPE with replaceable cartridges combined with automated injection to an LC/MS/MS system. Automated DBS extraction methods for a variety of analytes from different matrices are examined along with the use of different SPE cartridge sorbents. The resulting precision and accuracy data are provided.

and Determination of otic Autosampler and

ns were then introduced into the new Agilent C-MS/MS instrument.

UTION
y of sample handling steps are required prior to
sis of urine samples to accurately determine
concentrations. These steps typically begin
e enzymatic hydrolysis of analytes from their
sted forms to the native drug using enzymes
/beta-glucuronidase. The genetically modified,
eta-glucuronidase IMCzyme can hydrolyze
le drug classes within 30 minutes with high
ncy [2]. To ensure that the hydrolysis process
splete and reproducible, the pH, temperature
ration must be controlled and optimized for the
se used.

achieve the very low limits of detection necessary for compounds and their metabolites, it is often necessary to remove interfering matrix. Interferences are produced as a result of the hydrolysis procedure and naturally from the biological nature of the samples themselves. Solid phase extraction (SPE) is widely used, proven method for sample preparation and sample clean-up of hydrolyzed urine samples in the field of forensic analysis. Disposable Pipette Extraction (DPE) was developed as an alternative to traditional SPE, combining efficient and rapid extraction with significantly reduced solvent consumption. DPE relies

a reduction of the sample volume used from 1 to 0.5 mL serum and the use of smaller 1 mL format SPE cartridges. The analysis method has been validated according to GTFCh guidelines (Society of Toxicological and Forensic Chemistry). Limits of quantification below 1 ng/mL for THC and THC-OH, extraction efficiencies between 70 and 93 % and relative standard deviations between 3.3 and 10 % were achieved. The SPE system performs sample preparation in

Thank you for your attention!




1994

GERSTEL Inc.
USA



1967

GERSTEL GmbH & Co. KG
Germany



2004

GERSTEL K.K.
Japan



2001

GERSTEL AG
Switzerland

2010



GERSTEL LLP
Singapore



2018

GERSTEL Shanghai
Trading Ltd.
China

Thank you for your attention!

